



Forte offers the strong power that severe hearing loss users require and the sound quality that they appreciate. The Directional Microphone for conversations in noise environment and the multi-purpose Telecoil in this high performance, multi-memory product makes the Forte especially attractive to professionals, educators of the deaf and consumers.



- \* Volume Control
- \* Push button
- \* Directional Microphone
- \* 13A battery
- \* Telecoil
- \* Dimensions: 39.6x14.6x9.4 mm




# Forte

## Specifications

- \* High-Performance DSP hearing aid
- \* Up to 32 channels WDRC multi memories amplifier
- \* 12 Bands Frequency Shaping
- \* 12 Bands Layered Noise Reduction (Advanced)
- \* 3rd generation of industry-leading AFC technology
- \* Multi channels MPO Control
- \* Low Battery Warning
- \* Switch Tones
- \* Power-on Delay
- \* Adaptive Directional Microphone
- \* Up to 4 programs
- \* Programmable by CS44 Programming Cable


## Fitting Range



## Technical Data

	Forte
Maximum Output	136
HFA Avg Output	129
Maximum Gain	75
HFA Avg Gain	60
EQ Input Noise	≤32
THD@500Hz	≤3.0
THD@800Hz	≤2.5
THD@1600Hz	≤1.5
Freq Range F1	≤350
Freq Range F2	≥5000
Battery Current(mA)	≤2.0

## Testing Diagram



ANSI S3.22-2003

## IMHear Corporation

Add: 2711 Curtiss Street, Suite B, Downers Grove, IL, 60515, USA Tel: 1-630-395-9628  
 Website: [www.imhearing.com](http://www.imhearing.com) Email: [info@imhear.comcastbiz.net](mailto:info@imhear.comcastbiz.net)