




POCO

2-in-1 Hearing & Music Hearbuds



■ SELLING POINT

2-in-1 Hearing & Music ITE device

Indulge yourself with our comfortable earbud design to music streaming and hearing enhancement.

App self-tuning and remote-fitting

Connect wirelessly to iOS® or Android™ devices for SOUNDWEAR® APP remote control: Hearing check / Self tuning Sound Equalizer / Listening Mode Select

Long time use

The powerful battery capacity enable the POCO MU hearing aids to provide up to 8 hours of use. And the charging case will provide 4 times recharge cycle for both devices.

■ TECHNICAL DATA

| | |
|-------------------------------|----------|
| Maximum output OSPL 90(dBSPL) | 108 |
| HFA output OSPL 90(dBSPL) | 101 |
| Maximum Full-on gain(dB) | 35 |
| HFA Full-on gain(dB) | 28 |
| THD(%) | < 3 |
| EQ Input Noise(dB) | < 32 |
| Frequency Range (Hz) | 350-4500 |
| Battery Current (mA) | |

Data in accordance with ANSI S3.22-2003 IEC 60118-7:2005 2cc coupler. All technical(typical) data are subject to change without prior notice.

■ FEATURE

- EZCharge® : Magnetic contact charging
- Up to 32 channels WDRC
- Long use time
- Low battery warning
- Drop and Charge, Pick and On
- Intelligent noise reduction function
- User friendly operation button
- 3 listening programs
- In-ear detection sensor
- Resume of memory

■ PACKAGE INFO

- Type-C charging cable
- Silicone ear dome
- Charging case
- Brush

| | |
|-------------------------------|--------------|
| Battery type | Lithium |
| Working hours (Hr) Hearing | 8 hours |
| Working hours (Hr) Music | N/A |
| Charging time (per pair) | 4 times |
| Time of full charge (Devices) | 3 hours |
| Time of full charge (Base) | 2 hours |
| Size and weight (Device) | 20.8*19*16mm |
| Size and weight (Base) | 70*42.6*33mm |

IMHear Corporation

Add: 2711 Curtiss Street, Suite B, Downers Grove, IL, 60515, USA Tel: 1-630-395-9628
 Website: www.imhearing.com Email: info@imhear.comcastbiz.net