



# Poco X

2-in-1 Hearing & Music Hearbuds

		
TWS Earbud Design	Conversation enhancement	Music & Phone For All Platform



## SELLING POINT

### ALL PLATFORM CONNECTIVITY

Connect wirelessly to iOS® or Android™ devices both for music, phone and APP remote control.

### APP SELF-FIT

Connect wirelessly to iOS® or Android™ devices for SOUNDWEAR® APP remote control:  
Hearing Test / Self Adjustment  
Sound Equalizer / Listening Mode Select

### MUSIC STREAMING

Pamper yourself with our comfortable earbud design to music streaming and hearing enhancement.

## FEATURE

- Low Battery Consumption DSP
- Noise reduction
- Feedback cancellation
- In-Ear detection sensor
- Low battery warning
- Drop and Charge, Pick and On
- Popular earbud design for comfortable wearing
- Touch button for program switch (3 programs)
- Charging case with power bank
- Button and LEDs for battery charging status

## PACKAGE INFO

- Type-C charging cable
- Silicone ear dome
- Charging base

## TECHNICAL DATA

Maximum output OSPL 90(dBSPL)	117
HFA output OSPL 90(dBSPL)	112
Maximum Full-on gain(dB)	32
HFA Full-on gain(dB)	28
THD(%)	< 3
EQ Input Noise(dB)	< 28
Frequency Range (Hz)	350-5000
Battery Current (mA)	6

Battery type	Lithium
Working hours (Hr)   Hearing	6 hours
Working hours (Hr)   Music	6 hours
Charging time (per pair)	4 times
Time of full charge (Devices)	1.5 hours
Time of full charge (Base)	3 hours
Size and weight (Device)	30.5*23*23mm
Size and weight (Base)	62.0*59.6*30.3mm

Data in accordance with ANSI S3.22-2003 IEC 60118-7:2005 2cc coupler.  
All technical(typical) data are subject to change without prior notice.

## IMHear Corporation

Add: 2711 Curtiss Street, Suite B, Downers Grove, IL, 60515, USA Tel: 1-630-395-9628  
Website: www.imhearusa.com Email: info@imhear.comcastbiz.net