



Primo

Rechargeable Powerful Hearing Aids





Magnetic Charging Self-own tech Rechargeable



SELLING POINT

SMART CHARGING TECHNOLOGY

Our EZcharge® can prevent overheating during charging, reduce product power consumption, improve battery life and hearing aid safety.

LONG USE TIME

The powerful battery capacity enable the Primo P401 hearing aids to provide up to 55 hours of use. And the charging base offers 2 times recharge cycles for both hearing aids.

POWERFUL AND COMFORTABLE

Primo P401 offers maximum gain of 77dB which would fit profound hearing loss people, and also provides comfortable wearing.

FEATURE

- Up to 16 Channels WDRC SDP
- Super power BTE
- Adaptive Feedback Cancellation AFC Technology(Advanced)
- 10 Bands Layered Noise Reduction
- Wideband MPO Control
- Low battery warning
- Adaptive directional dual microphones
- 4 listening programs
- T-coil function

PACKAGE INFO

- Type-C charging cable
- Silicone ear dome
- Charging base
- Brush

TECHNICAL DATA

Maximum output OSPL 90(dBSPL)	135
HFA output OSPL 90(dBSPL)	129
Maximum Full-on gain(dB)	76
HFA Full-on gain(dB)	61
THD(%)	< 3
EQ Input Noise(dB)	< 32
Frequency Range (Hz)	450-3000

Battery Current (mA)	2.2 mA
Battery type	Lithium
Working hours (Hr) Hearing	55 hours
Time of full charge (Devices)	2-2.5 hours
Size and weight (Device)	9.1*18.3*34.8 mm
Size and weight (Base)	84.0*30*65.2mm

Data in accordance with ANSI S3.22-2003 IEC 60118-7:2005 2cc coupler. All technical(typical) data are subject to change without prior notice.

IMHear Corporation

Add: 2711 Curtiss Street, Suite B, Downers Grove, IL, 60515, USA Tel: 1-630-395-9628
Website: www.imhearing.com Email: info@imhear.comcastbiz.net