



Sax is an outstanding instrument that offers middle power for moderate hearing loss person. With stable and high performance, Sax offers a better choice for your cost-sensitive clients.




Sax

Specifications

- * Optimized DSP hearing aid
- * Up to 8 Channels WDRC Multi Memories Amplifier
- * 12 Bands Frequency Shaping
- * 10 Bands Layered Noise Reduction
- * Adaptive feedback cancellation
- * Up to four programs
- * Low Battery Warning


Fitting Range



Technical Data

	Sax
Maximum Output	125
HFA Avg Output	120
Maximum Gain	55
HFA Avg Gain	52
EQ Input Noise	≤30
THD@500Hz	≤3.0
THD@800Hz	≤3.0
THD@1600Hz	≤3.0
Freq Range F1	≤200
Freq Range F2	≥6000
Battery Current(mA)	≤0.8

Testing Diagram



ANSI S3.22-2003

IMHear Corporation

Add: 2711 Curtiss Street, Suite B, Downers Grove, IL, 60515, USA Tel: 1-630-395-9628
 Website: www.imhearusa.com Email: info@imhear.comcastbiz.net